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CYTOTOXIQUE	Indice de performance à la perméation
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	0
Cisplatine, 1.0 mg/ml (1,000 ppm)	5
Cyclophosphamide (Cytoxan), 20.0 mg/ml (20,000ppm)	5
Dacarbazine (DTIC), 10.0 mg/ml(10,000 ppm)	5
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	5
Etoposide (Toposar), 20.0 mg/ml (20,000 ppm)	5
Fluorouracil, 50.0 mg/ml (50,000 ppm)	5
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	5
Thiotepa, 10.0 mg/ml (10,000 ppm)	0
Methotrexate, 25 mg/ml, (25,000 ppm)	5
Mitomycin C, 0.5 mg/ml (500 ppm)	5
Vincristine Sulfate, 1.0 mg/ml (1,000 ppm)	5
Ifosfamide 50,0 mg/ml (50,000 ppm)	5
Mitoxantrone 2 mg/ml (2,000 ppm)	5

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DISPOSITIF MEDICAL

Testing. Development. Problem Solving.

March 24, 2011


• TEST REPORT •

PN 94809 A


CHEMICAL ANALYTICAL SERVICES

Prepared For:

Prepared By:

Tiffany L. Heller
Senior Technician

Approved By:

Ana C. Barbur, M.S.

Manager, Chemical Microbiological & Pharmaceutical Services

An A2LA Accredited Testing Laboratory — Certificate Numbers 255.01 & 255.02
ISO 9001:2008 Registered
Member of ACIL: The American Council of Independent Laboratories

ISO 9001:2008
Registered

Letters and reports are for the exclusive use of the clients to whom they are addressed and shall not be reproduced, except in full, without the written permission of Akron Rubber Development Laboratory, Inc. (ARDL). The information contained herein applies to the specific material, products or processes tested or evaluated. No warranty of any kind is herein constituted or implied. The liability of ARDL, Inc. shall be limited to the amount of consideration paid for services. ARDL, Inc. is accredited by A2LA for the test methods listed on the attached scope.

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DISPOSITIF MEDICAL

Page 1 of 4 – PN 94809 A

SUBJECT: Permeation testing per ASTM D 6978-05 on sample submitted by the above company.

RECEIVED: Glove sample identified as Non-Sterile Blue Powder Free Nitrile Examination Gloves – Cobalt Blue, Production Date June 10, 2010, Batch# 016105A27, Size Medium.

TESTING CHEMOTHERAPY DRUGS:

Table 1. List of the Testing Chemotherapy Drugs, Sources, and Expiration Dates

TESTING CHEMOTHERAPY DRUGS	CHEMICAL SOURCE
Carboplatin	Teva; Lot# 10D26LA; Expiration 04/2012
Carmustine (BCNU)	Bristol-Myers; Lot# 0E7004A; Expiration 05/2013
Cisplatin	Teva; Lot# 10G23KA; Expiration 01/2012
Cyclophosphamide (Cytoxan)	Sigma; Lot# 079K1569; Expiration 12/2011
Dacarbazine (DTIC)	Hospira; Lot# X022223AA; Expiration 05/2012
Doxorubicin Hydrochloride	Ben Venue; Lot# 1827762; Expiration 01/2012
Etoposide (Toposar)	Teva; Lot# 31311001B; Expiration 02/2013
Fluorouracil	APP; Lot# 6100345; Expiration 12/2011
Ifosfamide	Baxter; Lot# 0F337A; Expiration 06/2013
Methotrexate	Intas; Lot# K5340; Expiration 4/2011
Mitomycin C	Sigma; Lot# 048K1086; Expiration 1/2012
Mitoxantrone	Sigma; Lot# 050M1241; Expiration 12/2012
Paclitaxel (Taxol)	Hospira; Lot# W136865AB; Exp. 09/2011
Thiotepa	USP; Lot# I, Catalog# 66400; Expiration 12/2011
Vincristine Sulfate	USP; Lot# QOJ245; Expiration 08/2012

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Page 2 of 4 – PN 94809 A

COLLECTION MEDIA:

Table 2. Collection Media for Testing Chemotherapy Drugs

TEST CHEMICAL AND CONCENTRATION	COLLECTION MEDIUM
Carboplatin, 10 mg/ml (10,000 ppm)	Distilled Water
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	10% Ethanol Aqueous Solution
Cisplatin, 1.0 mg/ml (1,000 ppm)	Distilled Water
Cyclophosphamide (Cytoxan), 20.0 mg/ml (20,000ppm)	Distilled Water
Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm)	Distilled Water
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	Distilled Water
Etoposide (Toposar), 20.0 mg/ml (20,000 ppm)	Distilled Water
Fluorouracil, 50.0 mg/ml (50,000 ppm)	9.20 pH Sodium Hydroxide Solution
Ifosfamide, 50.0 mg/ml (50,000 ppm)	Distilled Water
Methotrexate, 25 mg/ml (25,000 ppm)	Distilled Water
Mitomycin C, 0.5 mg/ml (500 ppm)	Distilled Water
Mitoxantrone, 2 mg/ml (2,000ppm)	Distilled Water
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	30% Methanol Aqueous Solution
Thiotepa, 10.0 mg/ml (10,000 ppm)	Distilled Water
Vincristine Sulfate, 1.0 mg/ml (1,000 ppm)	Distilled Water

DETECTION METHOD OF CHEMICAL PERMEATION; UV/VIS ABSORPTION SPECTROMETRY:

Instrument: Perkin Elmer UV/VIS Spectrometer Lambda 25

UV/VIS Absorption Spectrometry was used to measure the absorbance of test chemicals which permeated through the specimens into the collection medium. The collection medium was circulated in a closed loop at 11 ml/minute of flow rate through the testing period. Data collection was performed according to the programmed schedule by means of UV Winlab software from the Perkin Elmer Corporation. The list of the characteristic wavelengths is shown below.

Table 3. Characteristic Wavelengths used in UV/VIS Absorption Spectrometry

TESTING CHEMOTHERAPY DRUGS	WAVELENGTH (nm)
Carboplatin	192
Carmustine	229
Cisplatin	199
Cyclophosphamide (Cytoxan)	200
Dacarbazine (DTIC)	320
Doxorubicin Hydrochloride	232
Etoposide (Toposar)	205
Fluorouracil	269
Ifosfamide (Ifex)	200
Methotrexate	303
Mitomycin C	217
Mitoxantrone	242
Paclitaxel (Taxol)	231
Thiotepa	199
Vincristine	220

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Page 3 of 4 – PN 94809 A

TESTING CONDITIONS:

Standard Test Method Used:	ASTM D 6978-05
Deviation From Standard Test Method:	Used 1" Permeation Cell
Analytical Method:	UV/VIS Spectrometry
Testing Temperature:	35.0°C ± 2.0
Collection System:	Closed Loop
Specimen Area Exposed:	5.067 cm ²
Selected Data Points:	25/test
Number of Specimens Tested:	3/test
Location Sampled From:	Cuff area
Comments/Other Conditions:	Magnetic stir bar was used in the sampling chamber

SAMPLE CHARACTERISTICS:

Table 4. Thickness characteristics for the tested specimens on: Non-Sterile Blue Powder Free Nitrile Examination Gloves – Cobalt Blue, Production Date June 10, 2010, Batch# 016105A27, Size Medium.

Testing Chemotherapy Drugs	Thickness (mm)			Average (mm)	Weight/Unit Area (g/m ²)
	#1	#2	#3		
Carboplatin	0.060	0.059	0.056	0.058	55.2
Carmustine	0.054	0.057	0.055	0.055	55.2
Cisplatin	0.056	0.055	0.055	0.055	55.2
Cyclophosphamide (Cytosan)	0.057	0.056	0.053	0.055	55.2
Dacarbazine (DTIC)	0.055	0.056	0.057	0.056	55.2
Doxorubicin Hydrochloride	0.052	0.056	0.056	0.055	55.2
Etoposide (Tosopar)	0.057	0.057	0.053	0.056	55.2
Fluorouracil	0.056	0.055	0.052	0.054	55.2
Ifosfamide (Ifex)	0.058	0.056	0.056	0.057	55.2
Methotrexate	0.056	0.053	0.056	0.055	55.2
Mitomycin C	0.054	0.055	0.057	0.055	55.2
Mitoxantrone	0.056	0.054	0.057	0.055	55.2
Paclitaxel (Taxol)	0.055	0.057	0.057	0.056	55.2
Thiotepa	0.052	0.059	0.060	0.057	55.2
Vincristine	0.057	0.056	0.052	0.055	55.2

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
Page 4 of 4 – PN 94809 A

RESULTS:

Table 5. Permeation Test Results on: Non-Sterile Blue Powder Free Nitrile Examination Gloves – Cobalt Blue. Production Date June 10, 2010, Batch# 016105A27, Size Medium.

TEST CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3) (Minutes)	AVERAGE STEADY STATE PERM. RATE (Specimen 1/2/3) ($\mu\text{g}/\text{cm}^2/\text{minute}$)	OTHER OBSERVATIONS
Carboplatin, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	1.82 (3.96, 1.82, 5.98)	2.4 (2.2, 2.4, 2.7)	Moderate swelling and no degradation
Cisplatin, 1.0 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Cyclophosphamide (Cytosan), 20.0 mg/ml (20,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Etoposide (Toposar), 20.0 mg/ml (20,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Fluorouracil, 50.0 mg/ml (50,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Ifosfamide, 50.0 mg/ml (50,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Methotrexate, 25 mg/ml (25,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Mitomycin C, 0.5 mg/ml (500 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Mitoxantrone, 2 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	No breakthrough up to 240 min.	N/A	Moderate swelling and no degradation
Thiotepa, 10.0 mg/ml (10,000 ppm)	0.93 (3.09, 5.38, 0.93)	1.9 (2.3, 2.2, 1.3)	Slight swelling and no degradation
Vincristine Sulfate, 1.0 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation

Tiffany L. Heller
Senior Technician
Pharmaceutical Services

AKRON RUBBER DEVELOPMENT LABORATORY, INC.

Ana C. Barbur, M.S.
Manager
Chemical, Microbiological & Pharmaceutical Services

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VIRUAL-1

Lowenkamp R&D Laboratories Ltd		Member:		ASTM Int'l., IEST, AZLA, ACS,	
ANALYTICAL RESEARCH				ADAC Int'l., ASOC, GAATW,	
P.O. Box 878, 1044 Lowenkamp Lane				IUPAC/IEU, NAIS	
Hazlehurst, MS 39083 U.S.A.				ASTM Lab #63450	
Tel/Fax: 601-894-2802 Mobile: 214-914-2276		Email: lionfight@aol.com			
LAB LOG NO. 12-18103 /2		CLIENT:		NO. 2	
DATE RECEIVED: June 29, 2012		CONTACT:			
DATE TESTED: July 3 to July 7, 2012		ADDRESS:		email: _____	
PHONE:		FAX: 05-679 2288		05-679 1188	
PRODUCT: Nitrile Exam Glove, On-Line Powder-Free		Tex. Finger		NBR (PF) F-T 8612F (C)	
SPECIFICATION: ASTM F 1671-07, VIRAL PENETRATION TEST - VIRAL BACTERIOPHAGE VIA Phi-X174		COLOR: Blue		SIZE: MEDIUM	
RECOVERY IS EXPECTED TO BE 100% +/- 2% ALLOWANCE		Broth Mix = 900 to 1200 PFU's Total Mix (40 to 50 PFU's/mL)		REPORTING LEVEL <1.0 PFU's/mL	
SAMPLE TO BE CUT 3" Dia. MINIMUM		PLAQUE FORMING UNITS (PFU) PER ML - >10 PFU/mL CONCENTRATION		BRAND:	
Batch/LOT No. NC216526B12		Pdn Date: 13/06/2012			
TEST: VIRUAL PENETRATION		0 kPa (0 psig) 5 minutes / 13.8 kPa (2 psig) 1 minute / 0 kPa (0 psig) 54 minutes, Procedure A			
GLOVES AGED AT:		24 HOURS @ 21 deg. C +/- 5 deg., Relative H 60%		Samples taken from palm of glove	
SAMPLE THICK & DESC.		SAMPLE THICK & DESC.		SAMPLE THICK & DESC.	
SEQ. BRAND NAME & DESC.		LOC. WT (g) NUMBER		LOT	
1 Nitrile Exam Glove, PF		PALM 0.080		SEE	
Blue		3" DIA 0.351		ABOVE	
2 Nitrile Exam Glove, PF		PALM 0.080		.	
Blue		3" DIA 0.385		.	
3 Nitrile Exam Glove, PF		PALM 0.080		.	
Blue		3" DIA 0.310		.	
7 CONTROL POSITIVE		sheet 0.145		.	
J45um nylon spheribach		3" DIA 0.288		.	
8 CONTROL NEGATIVE		sheet 0.072		.	
Mylar membrane		3" DIA 0.314		.	
I certify that the test results are performed to the required specification(s) and only reflect data obtained and/or observed from the samples provided for testing. The results do not reflect shipments prior to the stated PO or Lot Numbers and do not reflect the condition of future shipments. UNDER THE SPECIFICATION(S) APPLIED THE PRODUCT HAS BEEN FOUND TO BE:		ACCEPTABLE			
DATA REVIEWED BY:		NOTE: No problems observed.			
		Date: July 7, 2012			
William C. Lowenkamp Jr., Ph.D. Engr.		President			

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Signatory of EA, IAF and ILAC
Mutual Recognition Agreements

Sede operativa: 27029 VIGEVANO (PV) - C.so G. Brodolini, 18 - Tel. 0381.84722 - Fax 0381.73383 - E-mail: info@cimaonline.com - Internet: http://www.cimaonline.com

According to the EEC instruction 89/686 dated 21st of December 1989, concerning the standardisation of legislation of all Member Countries with regards to individual protection system and the relative legislative decree dated 4th of December 1992, N°. 475,

A.N.C.I. Servizi s.r.l. a socio unico - C.I.M.A.C. section
CENTRO ITALIANO MATERIALI DI APPLICAZIONE CALZATURIERA
Authorized with the decree issued by the Ministry of Industry of the Italian Republic
On the 11th of October 2000 – Community identification N°. 0465

grants:

CLOSURE REPORT OF EC CERTIFICATION

N° 0162/23853/16 - 00862C

MODULE C2 – EC TYPE CONFORMITY BASED ON INTERNAL PRODUCTION CONTROL WITH TESTS MADE UNDER CONTROL OFFICER AT RANDOM INTERVALS

For the following model of personal protective equipment of Cat. III:

**Protective glove against chemicals and micro organism - Article
"SENSISKIN bluple nitrile glove 240 mm"**

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The validity of the closure report shall be subject to surveillance activities provided for by Directive 89/686/EEC, to be implemented within 12 months from the issue date specified below.

Vigevano, 23 May 2017

Responsible of CE Certification of Glove
Dr. Sandro Milanesi

The Technical Responsible
Ing. Giuseppe Bellotti

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Signatory of EA, IAF and ILAC
Mutual Recognition Agreements

Sede operativa: 27029 VIGEVANO (PV) - C.so G. Brodolini, 19 - Tel. 0381.84722 - Fax 0381.73393 - E-mail: info@cimaonline.com - Internet: <http://www.cimaonline.com>

1. Description of personal protective equipment:

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Type of PPE: Protective glove against chemicals and micro organism

In compliance with EN 374: 2003 Part I II and III

Protection against anti neoplastic agents

Model of glove: glove with five fingers

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Signatory of EA, IAF and ILAC
Mutual Recognition Agreements

Sede operativa: 27029 VIGEVANO (PV) - C.so G. Brodolini, 19 - Tel. 0381.84722 - Fax 0381.73383 - E-mail: info@cimaonline.com - Internet: <http://www.cimaonline.com>

2. The tests and the examinations to verify the conformity of the article (in compliance with art. 11A of Directive 89/686/EEC – Decision 768/08/EC Module C2) are performed applying the following harmonized standards:

Standards applied:

- A - EN 420:2003 + A1: 2009 – General requirements for gloves.
- B - EN 374:2003 part I II e III – Glove against chemicals and micro organism
- C - EN 388:2016 – Protective gloves against mechanical risk.

3. The results of tests and examinations are contained in the following test reports:

C.I.M.A.C. 2017/0936 - RP - 2-RP-1 Dated 23 of May 2017

4. Requirements of the personal protective equipment:

Based on the test carried out, the model of PPE examined is conform to the model described in the CE type certificate and that there is uniformity in production activities.

The model of Protective glove against chemicals and micro organism - Article "SENSISKIN bluple nitrile glove 240 mm" conforms:

- the requirements of EN 420:2003 + A1 2009 standard, points 4.1, 4.2, 4.4, 5 and 5.2;
- the following permeation levels as specified in EN 374:2003 Part I standard;

	Performance level
Diethyl amine	2
Sodium Hydroxide solution 40%	2
Sulphuric acid solution 96%	2

- the performance levels as specified in prospectus 1 of EN 388:2016 standard;

	Performance level
6.1 Abrasion resistance	0
6.2 Cut resistance	0
6.3 Tear resistance	0
6.4 Puncture resistance	0
5.2 Dexterity	5

In the model of Personal Protective Equipment and its components was not detected the presence of dangerous substances listed in Annex XVII of Regulation 1907/2006/EC and subsequent amendments and additions.

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Signatory of EA, IAF and ILAC
Mutual Recognition Agreements

Sede operativa: 27029 VIGEVANO (PV) - C.so G. Brodolini, 19 - Tel. 0381.84722 - Fax 0381.73393 - E-mail: info@cimaonline.com - Internet: <http://www.cimaonline.com>

5. Marking of the personal protective equipment:

The following information is provided on the box containing the gloves:

- the "CE" mark
- the identification of Competent Body: 0465
- the article code: SENSISKIN bluple nitrile glove 240 mm
- the manufacturer name: LABORATOIRES EUROMEDIS
- the glove size
- the symbols concerning the protection provided: Protective gloves against chemicals and micro-organisms.
- the permeation performance levels as specified in EN 374:2003 Part I standard:

	Performance level	Code
Diethyl amine	2	G
Sodium Hydroxide solution 40%	2	K
Sulphuric acid solution 96%	2	L

- the performance levels as specified in prospectus 1 of EN 388:2016 standard;

	Performance level
6.1 Abrasion resistance	0
6.2 Cut resistance	0
6.3 Tear resistance	0
6.4 Puncture resistance	0
5.2 Dexterity	5

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Signatory of EA, IAF and ILAC
Mutual Recognition Agreements

Sede operativa: 27029 VIGEVANO (PV) - C.so G. Brodolini, 19 - Tel. 0381.84722 - Fax 0381.73393 - E-mail: info@cimaonline.com - Internet: <http://www.cimaonline.com>

6. Notes:

- "*making available on the market*" shall mean any supply of a product for distribution, consumption or use on the Community market in the course of a commercial activity, whether in return for payment or free of charge.
- "*placing on the market*" shall mean the first making available of a product on the Community market.
- "*manufacturer*" shall mean any natural or legal person who manufactures a product or has a product designed or manufactured, and markets that product under his name or trademark.
- "*authorised representative*" shall mean any natural or legal person established within the Community who has received a written mandate from a manufacturer to act on his behalf in relation to specified tasks.
- "*harmonised standard*" shall mean a standard adopted by one of the European standardisation bodies listed in Annex I to Directive 98/34/EC on the basis of a request made by the Commission in accordance with Article 6 of that Directive.
- "*accreditation*" shall mean an attestation by a national accreditation body that a conformity assessment body meets the requirements set by harmonised standards and, where applicable, any additional requirements including those set out in relevant sectoral schemes, to carry out a specific conformity assessment activity.
- "*conformity assessment*" shall mean the process demonstrating whether specified requirements relating to a product, process, service, system, person or body have been fulfilled.
- "*recall*" shall mean any measure aimed at achieving the return of a product that has already been made available to the end user.
- "*withdrawal*" shall mean any measure aimed at preventing a product in the supply chain from being made available on the market.
- "**CE** marking" shall mean a marking by which the manufacturer indicates that the product is in conformity with the applicable requirements set out in Community harmonisation legislation providing for its affixing.
- The content of this EC Type-Examination Certificate is referred to the tested personal protective equipment only.
- This EC Type-Examination Certificate may be integrally duplicated; the copy must be faithful, legible (if pint size) and must contain the bold caption "TRUE COPY".

2017/0936-2-CEGC-2

Closure report of EC Certification N° 0162/23853/16 - 00862C dated 23/05/2017 - Page: 5 of 5 – M52 Rev.0 06.03.17

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SENSISKIN © Gant nitrile sans poudre Bleu-Violet