



CHEMICAL RESISTANCE

A NITRITEX INFORMATION GUIDE



GB10/81851



GB04/61092

UK e-mail: info@nitritex.com
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Nitritex Ltd, CB8 7YY, United Kingdom

Disposable Gloves Chemical Resistance Guide

Breakthrough time in minutes

<15 minutes

>15 but <180 minutes

>180 minutes

Chemical	Nitrile	Neoprene	Latex	Vinyl
Acetaldehyde	<0.6	10.2	7	1
Acetic Acid	5	360	21	45
Acetic Anhydride		210	3	
Acetone	3	2.4	2.4	<1
Acetonitrile (Methyl Cyanide)	<5	<10.8	<0.6	<1
Acrolein	4.2			
Acrylic Acid		70	80	60
Acrylonitrile	<5			
Aldehyde	4			
Allyl Alcohol		94.8		60
Allylamine			<1.2	
Ammonium Fluoride	>360	>360	>360	>360
Ammonium Hydroxide	360	360	90	240
Amyl Acetate (Isoamyl Acetate)	<5	5.4	5.4	<1
Amyl Alcohol	30	321	7.2	12
Amyl Nitrile	175.8	46.8		
Aniline	18	30	25	62
Benzaldehyde	<5	39	10	<1
Benzene	4.2	1.2	0.6	<1
Benzenesulfonic Acid		>1200		
Benzethonium Chloride		>480	>480	
Benzonitrile (Phenyl Cyanide)			<0.6	
Benzoyl Chloride		15		
Bisphthalate	259.8	120	>360	
Boric Acid	>480	>480		
Bromobenzene	13.2			
Bromoethanol			1.2	
Bromopropanol		>480	1.2	
Bromopropionic Acid		180	190	180
Butadiene		46.8	<1.2	
Butanol (Isobutanol)	30	10	1.2	
Butyl Acetate	<5	3.6	1.8	<1
Butyl Acrylate	67.8			
Butylamine	19.8	12	1.2	
Butyl Cellosolve (Butoxyethanol)	9	90	45	
Butyl Chloride (Chlorobutane)	12			
Butyl Nitrite	97.8			
Butyltoluene	>360	73.2		
Butyraldehyde		25.2		
Butyrolactone		10	60	
Carbon Disulfide	1		<5	<1

Chemical	Nitrile	Neoprene	Latex	Vinyl
Carbon Tetrachloride	5	4.8	3.6	25
Cellosolve (Ethoxyethanol)		45	25	
Cellosolve Acetate (Ethoxyethyl Acetate)	<5	25	10	
Chlorine		>480		
Chloroaceto-nitrile			<0.6	
Chlorobenzene	<5	10.8	<5	<5
Chlorobutadiene (Chlorophrene)	3.6	3		
Chloroethanol (Ethylene Chlorohydrin)		298.8		
Chloroform	2.4	0.6	0.6	<1
Chloronaphthalene	174		<5	<5
Chloronitropropane			1.2	
Chloropropanol			<0.6	
Chlorotoluene (Benzyl Chloride)	15			
Chromic Acid	240	75	70.2	
Citric Acid		>360	>360	>360
Copper		>360		>360
Creosote		270		
Cresol	<5	>60	13.8	
Crotonaldehyde		21		
Cyclohexane	360	6	1.8	
Cyclohexanol	360	150	10	360
Cyclohexanone	<5			<1
Cyclohexylamine		36	1.2	
Cyclopentanone				
Decanal (Decyl Aldehyde)		240		
Diacetone Alcohol (Hydroxy Methyl Pentanone)		300	15	<1
Diamylamine	>480	129		
Dibromoethane (Ethylene Dibromide)	27	4.8	<1.2	
Dibutylamine	>480			
Dibutyl Phthalate (Butyl Phthalate)	30	120	16.8	<1
Dichlorobenzene	<5			
Dichlorobutene	2.4	10.2		
Dichloroethane (Ethylene Dichloride)	2.4	1.8	0.6	
Dichloroethylene	7.2			
Diesel	>240			
Diethanolamine	>480	>480		
Diethylamine	<5		<5	
Diethylamino-ethanol	>480			
Diethylenedioxiide (Dioxane)	<5	6.4	2.4	
Diethylenetriamine		>480		
Diethylphthalate				
Diisobutylamine	>480	52.2		
Diisobutyl Ketone	120	15	15	

Chemical	Nitrile	Neoprene	Latex	Vinyl
Diisopropylamine	195	40.2		
Dimethylacetamide	<5		15	<5
Dimethylamine		>480	1.8	
Dimethylamino-propylamine		28.8	0.6	
Dimethylbutylamine	81			
Dimethyl-ethanolamine (Dimethyl-aminoethanol)	>480	235.2	4.8	
Dimethylformamide	<5	1.2	25	<1
Dimethylhydrazine	6	37.8	3.6	
Dimethylmercury			0.25	
Dimethylsulfoxide	28.2	60	79.8	<15
Dimethyl-vinylchloride	9			
Diocylphthalate		>360	<5	<5
Divinylbenzene	60			
Dodecane				
Epichlorohydrin	19.8	15	<1.2	<1
Epoxybutane (Butylene Oxide)		4.2		
Ethanol	24	49.2	12	
Ethanolamine	360	360	210	
Ether (Ethyl Ether Diethyl Ether)	13.8	10	10.2	
Ethidium Bromide				
Ethyl Acetate	<5	12	4.8;	<1
Ethyl Acrylate		48		
Ethylamine	66			
Ethylbenzene	<5			
Ethylbromide		4.2		
Ethylbutylamine		73.2		
Ethyl Cellosolve (Ethoxyethanol)	91.8	244.8	<0.6	
Ethyl Cyanide			<0.6	
Ethylenediamine		399	4.8	
Ethylene Glycol	360	360	360	>360
Ethyleneimine		<4.8		
Ethylglycol Ether		45	25	<15
Ethylhexanoic Acid	>240	>240		
Ethylhexanol		>480		
Ethyl Methacrylate	22.8			
Formaldehyde	>1260	120	6	80
Formic Acid	5	>360	120	360
Freon (Dichlorodifluoro-methane)	10.2	3	2.4	2
Furaldehyde (Furfural)	<5	19.8	15	<5
Gasoline	30		<5	<5
Glutaraldehyde	>240	>480		>360
Heptane	360	45	1.2	
Hexachloro-cyclopentadiene	>480			
Hexamethyl-disilazane		50	15	<15
Hexamethyl-phosphoramidate	90			

Chemical	Nitrile	Neoprene	Latex	Vinyl
Hexane	78.6	3.6	4.8	<2
Hydraulic Fluid	>240		<5	
Hydrazine	>480	>960	150	>360
Hydrochloric Acid (Muriatic Acid)	360	360	290	300
Hydrofluoric Acid	120	60	90	
Hydrogen Peroxide	>360	4.8	>480	>360
Hydrogen Phosphide		10.2	30	
Hydroquinone	>360	>360	>360	>360
Imino-bispropylamine		>480	6	
Isoprene	52.2	16.2		
Kerosene	>360	>360	<5	360
Lactic Acid	>360	>360	>360	>360
Lauric Acid	>360	>360	>360	15
Limonene		64.8		
Maleic Acid	>360	>360	>360	>360
Methacrylic Acid	10.2			
Methacrylonitrile	7		<1.2	
Methanesulfonic Acid		>240		
Methanol	10.8	15	1.8	45
Methoxymethyl-pentanone		99		
Methoxypropanol				
Methyl acetate			<1.2	
Methyl Acrylate		15	1.2	
Methylamine	>480	270	25.2	
Methylamino-propylamine		63	3	
Methyl Bromide			<5	
Methyl Butyl Ether	5		<5	<5
Methyl Cellosolve	40.2	25	20	
Methyl Chloride			0.6	
Methylene Dichloride (Methylene Chloride)	1.8	0.6		<0.5
Methylethanolamine		>480		
Methyl Ethyl Ketone	3.6	2.4	1.2	<1
Methyl Ethyl Ketone Peroxide		>240	45	
Methyl Glycol Ether		25	20	
Methylhexanone (Methyl Isoamyl Ketone)	<5			
Methyl Iodide	0.6	0.6	1.8	<1
Methyl Isobutyl Ketone	12	15	6	<5
Methyl Isocyanate		0.6	0.6	
Methyl Methacrylate	<5		<1.2	<1
Methylpyrrolidone	<5		75	
Mineral Oil	>240			150
Mineral Spirits	>360	90	<5	150
Monoethanolamine	360	360	50	>360
Monoethylamine	66			
Morpholine	<5		20	<5
Naphtha	>360	15	<5	120
Nickel		>360		

Chemical	Nitrile	Neoprene	Latex	Vinyl
Nitric Acid	5	79.8	<5	110
Nitrobenzene	<5	40.2	4.8	<5
Nitroethane		49.2	1.8	1
Nitrohydrochloric Acid (Aqua Regia)		45	<5	
Nitromethane	30	60	<1.2	<1
Nitropropane	12	5	1.8	<1
Nonylphenol		>1200		
Octane (Isoctane)	360	60	<5	
Octanol	360	360	30	
Oleic Acid	>360	60	30	90
Oxalic Acid	360	360	360	>360
Palmitic Acid	30	>360	4.8	75
Pentachlorophenol	>780	6	<5	180
Pentane	1.8	6.6	0.6	<1
Perchloric Acid	360	360	360	>360
Perchloroethylene			<5	<1
Petroleum Ether	>240			
Phenol	31.8	40.2	16.2	75
Phenolphthalein	>480	>480	>480	
Phosphoric Acid	360	360	360	>360
Phosphorus Oxychloride		<0.6		
Picric Acid		150	<5	40
Polychlorinated Biphenyl		1440	4.8	
Potassium Hydroxide	360	180	79.8	>360
Promethazine Hydrochloride	>480			
Propanol (Isopropanol)	30	90	7.2	90
Propanolamine (Monoisopropanolamine)		>480	30	
Propiolactone			19.8	
Propionaldehyde		12		
Propylacetate	16.8		4.8	<2
Propylamine		13.8		
Propylenediamine (Diaminopropane)		271.9	3	
Propyleneglycol			>360	
Propylene Oxide			<0.6	<1
Propyl Ether (Isopropyl Ether)	>60	42.6	3.6	
Propylmethacrylate	60			
Pyridine	5.4	1.8	2.4	<1
Rubber Solvent		30	<5	<5
Silicon Etch		>360	<5	150
Sodium Hydroxide	360	360	360	>360
Sodium Hypochlorite	360	360	360	
Stoddard Solvent	>240	180	<5	57
Styrene	30	12	10.2	10
Sulfuric Acid (Oleum)	5	70.2	<5	26
Tannic Acid	>360	>360	>360	>360
Tetrachloroethane	13.2	5.4	1.8	
Tetrachloroethylene	5	6	<1.2	

Chemical	Nitrile	Neoprene	Latex	Vinyl
Tetraethylene-pentamine		>480	106	
Tetrafluoroethylene		>480		
Tetrahydrofuran	0.6	1.2	1.2	<0.5
Tetramethylenediamine	108			
Toluene	<5	1.2	0.6	<0.5
Toluene Diisocyanate	222		7	5
Toluenesulfonic Acid		>480		
Triallylamine	>480	63		
Trichloroacetonitrile		67.2		
Trichlorobenzene	<5	60	4.8	
Trichloroethane	1.8	2.4	1.2	
Trichloroethylene	<5	1.8	0.6	<1
Trichloropropane	21			
Tricresylphosphate	60	>360	45	>360
Triethanolamine	>480	>360	60	>360
Triethylamine	>480	37.2		
Triethylene Tetramine	>480	>480		
Trifluoroethanol	7.2	>60	>60	
Tripopylamine	>480	>480		
Turpentine	30		<5	<5
Valeronitrile		40.8	1.8	
Vinyl Chloride	342			
Vinylcyclohexane	391.8			
Vinylidene fluoride		<1.2	<1.2	
Xylene	<5	3	1.2	<1

Note: The permeation rates are intended only as a guide. The suitability of all gloves must be determined by the user through their own testing.

This guide should not be considered as a warranty that any glove is fit for purpose.

www.nitritex.com

For more information please contact your nearest Nitritex office.



GB04/61092



GB10/81851

UNITED KINGDOM (Headquarters) info@nitritex.com
MALAYSIA (Manufacturing) info@nitritex.com.my
CANADA info@nitritexcanada.com
SINGAPORE info@nitritex.com.sg

