

# Contec Filtered Alcohols

70% IPA or Denatured Ethanol blended with purified water

Contec alcohols are a blend of 70% v/v Isopropanol or denatured ethanol with 30% purified water. The alcohol blend is 0.2 micron filtered, filled and bagged in a Grade C (ISO Class 7) cleanroom.

Contec's ethanol is only denatured with IPA creating a very low residue product.

Supplied as either trigger sprays fitted with a protected system, which ensures sterility throughout use, or 5L capped containers for when larger volumes are required.



| FEATURES   | BENEFITS   |
|--|--|
| GMP manufactured under Grade A air flow in a Grade C cleanroom | • Ensures the alcohol, container and packaging are free from contamination and particulates  |
| No denaturants other than IPA                                  | • Completely residue free  |
| Trigger sprayer and "bag in bottle" protected system           | • Prevents bottle contents from becoming contaminated during use<br>• No need to discard unused product during shelf life<br>• Bottles can be completely emptied eliminating wastage |
| Trigger sprayer can be set to jet or spray                     | • Large droplet size reduces the risk of inhalation and provides good surface coverage   |
| Double bagged packed in linear tear packaging                  | • Each bag is easy to open even when wearing gloves<br>• Facilitates transfer disinfection into cleanroom  |

## Product Specification

- 70% v/v IPA or denatured ethanol with 30% purified water
- Ready to use with a 3 year shelf life
- 0.2 micron filtered
- Double bagged in polyethylene linear tear bags
- Tamper evident caps and triggers
- Certificates of analysis provided with every batch.
- Full supporting documentation including efficacy information (EN1276 and EN1650) is available on request.



For more information or to request a sample, email [infoeu@contecinc.com](mailto:infoeu@contecinc.com) or phone + 33 (0) 2 97 43 76 98.

Copyright © 2019 Contec, Inc. All rights reserved. INT119 110419

## Product Application

Ready to use so can be applied directly to the surface ensuring even coverage. Physically remove contaminants from the surface using an appropriate cleanroom wipe and a recommended wiping technique for optimum contamination control.

| Part No. | Description  | Size | Packaging |
|----------|--|------|-----------|
| FBT170I  | Contec IPA<br>70% IPA with purified water, trigger spray                           | 1L   | 6 x 1L    |
| FBC570I  | Contec IPA<br>70% IPA with purified water, capped                                  | 5L   | 2 x 5L    |
| FBT170DE | Contec Denatured Ethanol<br>70% denatured ethanol in purified water, trigger spray | 1L   | 6 x 1L    |
| FBC570DE | Contec Denatured Ethanol<br>70% denatured ethanol in purified water, capped        | 5L   | 2 x 5L    |

Contec offers a range of innovative wipes and mopping systems designed for efficient and effective cleaning in controlled environments. See [www.contecinc.com](http://www.contecinc.com) or contact your local sales representative for advice and information about the most suitable products and cleaning techniques for your facility.

*Use biocides safely. Always read the label and product information before use.*